

Stop Wasting Milk Because of Staph. Aureus Mastitis





The ProScience ProStaph I™ mastitis test is the newest management tool available for milk producers concerned about detecting and controlling Staph. aureus.

Why use the ProScience ProStaph I[™] test? Convenience and accuracy. ProStaph I[™] uses either preserved or raw milk samples. The test is highly accurate and you avoid the hassle of swabbing teats for sterile culture samples. There is no interference with your milking routine.

Use ProStaph I[™] and receive accurate management information you and your veterinarian need to help you manage for top production and maximum milk quality premiums.

Call Pennsylvania DHIA (800) DHI-TEST or ProScience Corporation (800) 658-8868